



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Power Retail Prescription Form

JANUARY 2025
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX40 HD with Modulite HD seating

Standard Features

- Max user weight 160kg (25 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite HD adjustable seating system
- HD Telescopic Seat Frame and Plate with Side Rails width 530 to 550 mm with an additional 50mm, (2") by arms
- Seat depths 510mm - 580mm stepless adjustment
- Tension adjustable backrest
- Standard T-Armrests, detachable with short armpad
- Manual settable seat Tilt (0° - 20°)
- Manual settable backrest Recline (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX control REM110 with standard bracket
- Choice of 5 Shroud colours, 5 Rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension Control - Traction - Comfort
- 6 km/h (4mph) 4 pole and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £6,234



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
JBH0002	AVIVA RX 40 HD Base Model Max user weight up to 160 kg HD Version, includes RX 40 Performance Package Wide Base only	£6,234
JBX1300	Special Extended lead time: + 10 days	POA
CHASSIS WIDTH*		
JBH2830	Chassis Wide (635 mm/25")	◆

SEATING TYPE*			
JBX2101	<input type="checkbox"/>	Invacare® Modulite HD Seat Unit On RX 40 HD only Tension Adjustable Upholstry only JBX2171	(NC 2273) <input type="radio"/>
SEAT FRAME*			
JBX2102	<input type="checkbox"/>	HD Telescopic Seat Frame and Plate with Side Rails width 530 to 550 mm	<input type="radio"/>
SEAT WIDTH*			
JBX2124	<input type="checkbox"/>	530 mm Seat Plate	<input type="radio"/>
JBX2104	<input type="checkbox"/>	550 mm Seat Plate	<input type="radio"/>
SEAT DEPTH*			
JBX2136	<input type="checkbox"/>	Seat Depth 510 mm	<input type="radio"/>
JBX2137	<input type="checkbox"/>	580 mm	<input type="radio"/>
SEAT CUSHION*			
JBX2148	<input type="checkbox"/>	Seat Cushion (90 mm high)	<input type="radio"/>
JBX2140	<input type="checkbox"/>	Without seat cushion	<input type="radio"/>
BACKREST SUPPORT*			
JBX2162	<input type="checkbox"/>	Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrix Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	<input checked="" type="radio"/>
BACKREST UPHOLSTERY / CUSHION*			
JBX2171	<input type="checkbox"/>	Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm MUST BE SELECTED WITH HD SEATING JBX2102	(NC 2230, 2247) <input checked="" type="radio"/>
UPHOLSTERY TYPE*			
JBX2156	<input type="checkbox"/>	Black Nylon	<input checked="" type="radio"/>
ARMREST (Pair)*			
JBX2191	<input type="checkbox"/>	T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	<input type="radio"/>
JBX2181	<input type="checkbox"/>	T-Arm detachable, Extended height Adjustment range from seat plate: 295 - 365 mm	£41
ARMPAD (Pair)*			
JBX2195	<input type="checkbox"/>	Short Armpad (325 mm)	<input checked="" type="radio"/>
JBX2196	<input type="checkbox"/>	Long Armpad (375 mm)	<input type="radio"/>
JBX2197	<input type="checkbox"/>	Long Ergonomic Armpad (397 mm)	(NC 2230, 2247) <input type="radio"/>
HEADREST / NECKREST - Adequate headrest option must be selected*			
JBX1350	<input type="checkbox"/>	Neckrest depth adjustable	(NC 2230, 2247) £81
JBX1432	<input type="checkbox"/>	Headrest sideways adjustable	(NC 2230, 2247) £226
JBX1427	<input type="checkbox"/>	Headrest with cheek support	(NC 2230, 2247) £226
JBX2324	<input type="checkbox"/>	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware	(NC 2230, 2247) £361
JBX2325	<input type="checkbox"/>	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23) Complete with Elan Headrest Hardware	(NC 2230, 2247) £361
JBX2327	<input type="checkbox"/>	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard Complete with Elan Headrest Hardware	(NC 2230, 2247) £489
JBX1495	<input type="checkbox"/>	Non Invacare Headrest - no mounting hardware With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	<input type="radio"/>

SEAT LIFTER / SEAT TILT*				
JBX1091		Manually settable Seat Tilt 0 - 15°	(NC 1401)	○
JBX1092		Powered Tilt (0° - 20°)	(NC 2201)	£404
JBX1093		Powered Tilt (0° - 30°) with CoG shift	(NC 2201)	£1,073
JBX1094		Seat lifter with powered tilt (300 mm, 30°) with CoG shift, 6 point Tie downs Must be selected	(NC 2201, 2220)	£1,893
SEAT TO GROUND HEIGHT*				
JBX2857		Seat-to-Floor Height 435 mm		○
JBX2858		Seat-to-Floor Height 460 mm		○
JBX2859		Seat-to-Floor Height 485 mm		○
BACKREST RECLINER*				
JBX0975		Manual Settable Recliner (in 3.75° steps) 90° - 120°		◆
JBX0977		Powered Recliner 30°	(NC 2201)	£505
CENTRE MOUNTED LEGREST*				
JBH0556		Centre mount legrest, angle settable Length: 280 - 385 mm	(NC 0517, 1244, 1416)	£373
JBH0589		LNx Power centre mount legrest, length compensating Length: 340 - 410 mm	(NC 2201, 1244, 1416)	£1,640
JBH2855		Without Legrests		○
CENTRE MOUNT FOOTPLATES / FOOT PLATFORM*				
JBH2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")		£103
JBH0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")		○
JBH0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5")		£135
JBH0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")		£270
JBX3001		Please select If Centre mount footplate is not required.	(NC 0589)	○
FOOTREST HANGER TYPE*				
JBH0503		Hanger Vari F pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with tools		◆
JBH0533		Hanger Vari A pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with handwheel and calf pad		£233
JBH0531		Hanger ADE pair Lower leg length: 290 - 460 mm Powered Angle adjustment: 80 - 0° with calf pad	(NC 2201)	£669
JBH0501		Without Hanger pair	(NC 0504, 0513)	○
FOOTPLATE TYPE*				
JBH0512		Fixed Footplate pair	(NC 0533, 0531)	◆
JBH0504		Angle Adjustable Footplate pair		£32
JBH0513		Angle and depth adjustable Footplate pair		£118
JBX3003		No Side Mounted Footplate required	(NC 0503, 0533, 0531)	○
CONTROLS*				
JBH2201		Remote LiNX DLX-REM110 = Drive only 1 Drive function with Drive only power module	(NC 2837, 1092, 1093, 1094, 0977, 0589, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2258, 2219, 2245, 2270, 1908)	◆
JBH2203		Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators	(NC 1437, 2231, 2232, 2233, 2230, 2247, 2257, 2240, 2270, 2258, 1908, 1909)	○
JBH2204		Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights	(NC 2232, 2233, 2230, 2247, 2257, 2240, 2270, 2258)	£151
JBH2205		Rehab Remote LiNX, touch-screen user interface DLX-REM400 Configurable jack sockets	(NC 2220)	£341
JBH2206		Rehab Remote LiNX without joystick for alternative access users DLX-REM500 Requires a user input/alternative joystick: 2231, 2232, 2229, 2230, 2247, 2233	(NC 2220, 2270, 1207, 1208, 1707, 1708, 1709, 1710)	£776

REM400 TOGGLE OPTION			
JBH2270		Assignable Toggle Switch Module f. REM400 (DLX-SW400TGL)	(NC 2203, 2204, 2206) £316
CONTROL BRACKETRY*			
JBH1207		Standard Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404) ○
JBH1208		Standard Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405) ○
JBH1707		Swing Away Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404) £133
JBH1708		Swing Away Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405) £133
JBH1709		Matrx Resolve Swing in/out bracket RH	(NC 2283, 2287, 2207, 2265, 2275, 1404) £219
JBH1710		Matrx Resolve Swing in/out bracket LH	(NC 2284, 2288, 2208, 2266, 2276, 1405) £219
JBH2213		Centre Mount, swing-away to RH, Nucleus	(NC 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404, 1405) £498
JBH2214		Centre Mount, swing-away to LH, Nucleus	(NC 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1404, 1405) £498
JBH2227		Bracket RH, mounted for Attendant	(NC 2220, 2128, 2219, 1221, 1312) £64
JBH2228		Bracket LH, mounted for Attendant	(NC 2220, 2128, 2219, 1222, 1363) £64
JBX2239		Individual mounting option per request (Special) Requires additional information/specification	POA
DUAL CONTROL			
JBH2220		Dual Control LiNX (DLX-REM050 as secondary remote)	(NC 2273, 2205, 2206, 1221, 1222, 2227, 2228) £395
JBH2219		Compact Dual Control LiNX (DLX-ACU200)	(NC 1221, 1222, 2227, 2228) £395
INVACARE LiNX IDC DUAL CONTROL			
JBH2245		Invacare LiNX IDC Dual control	(NC 2259, 2260) £750
DUAL CONTROL BRACKETRY			
JBH2259		Bracket RH, Attendant	(NC 2291, 2245, 1312) ○
JBH2260		Bracket LH, Attendant	(NC 2292, 2245, 1363) ○
INVACARE LiNX IDC DUAL CONTROL BRACKETRY			
JBH1221		Bracket RH, detachable	(NC 2291, 2220, 2219, 2227, 1312) ○
JBH1222		Bracket LH, detachable	(NC 2292, 2220, 2219, 2228, 1363) ○
LiNX COMPACT REMOTE			
JBH2231		Proportional Compact Remote DLX-CR400	(NC 2203, 2204) £719
JBH2232		Proportional Compact Remote Low Force DLX-CR-400LF	(NC 2203, 2204) £960
LiNX COMPACT REMOTE BRACKET			
JBH2207		RH Mount	(NC 1707, 2213, 2217, 1404) ○
JBH2208		LH Mount	(NC 1708, 2214, 2218, 1405) ○
JBH2265		Centre Mount, swing-away to RH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!	(NC 1707, 2213, 2217, 1404, 1405, 1415) £498
JBH2266		Centre Mount, swing-away to LH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!	(NC 1708, 2214, 2218, 1404, 1405, 1339) £498
JBX2215		Individual mounting option per request (Special) Requires additional information/specification	POA
CHIN CONTROL			
UNI9250		Mo-vis Micro joystick (Please specify in Special Comments box mounting request if required)	POA
UNI9251		Mo-vis Multi joystick (Please specify in Special Comments box mounting request if required)	POA
LiNX HEAD CONTROL			
JBH2230		Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M3	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2171, 2197, 1350, 1432, 1427) £3,719
JBH2247		Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2171, 2197, 1350, 1432, 1427) £4,356

LINX FOUR-SWITCH PROXIMITY ARRAY			
JBH2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2230, 2247, 2233, 2240, 2128) £5,760
JBH2275		Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1707, 2213, 2217, 1404, 1405, 1415) £498
JBH2276		Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1708, 2214, 2218, 1404, 1405, 1339) £498
JBX2269		Individual mounting option per request (Special) Requires additional information/specification	POA
LINX WIRELESS ACCESSORIES FOR HEAD CONTROL AND PROXIMITY CONTROL			
JBH2249		Wireless Triple Switch Receiver (ASL557-3)	£341
JBH2250		Wireless Mouse Emulator (ASL558)	£817
LINX MICRO EXTREMITY CONTROL BRACKET			
JBH2233		MEC Proportional Joystick (Micro Extremity Control ASL130)	(NC 2203, 2204, 2230, 2247, 2257, 2128) £3,612
LINX HEAVY DUTY REMOTE			
JBH2240		Proportional Heavy Duty Remote Joystick w mode switch (ASL133)	(NC 2203, 2204, 2230, 2247, 2257, 2128) £2,013
LINX MICRO EXTREMITY AND HEAVY DUTY CONTROL BRACKET			
JBH2285		Tray Mount RH (100x150 mm)	(NC 1707, 2213, 2217, 1404, 1405, 1415) ○
JBH2286		Tray Mount LH (100x150 mm)	(NC 1708, 2214, 2218, 1404, 1405, 1339) ○
JBH2283		Centre Mount, swing-away to RH, Nucleus	(NC 1707, 2213, 2217, 1404, 1405, 1415) £498
JBH2284		Centre Mount, swing-away to LH, Nucleus	(NC 1708, 2214, 2218, 1404, 1405, 1339) £498
JBH2287		Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1707, 2213, 2217, 1404, 1405, 1415) £498
JBH2288		Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1708, 2214, 2218, 1404, 1405, 1339) £498
JBX2289		Individual mounting option per request (Special) Requires additional information/specification	POA
LINX 10-WAY SWITCH			
JBH2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access	£505
LINX 10-WAY SWITCH BRACKET			
JBH2217		Bracket RH, mounted for user	(NC 2207, 2287, 2283, 2275, 2213, 2265, 1404) ○
JBH2218		Bracket LH, mounted for user	(NC 2208, 2288, 2284, 2286, 2276, 2214, 2266, 1405) ○
JBH2291		Bracket RH, mounted for attendant	(NC 2259, 1221) £64
JBH2292		Bracket LH, mounted for attendant	(NC 2260, 1222) £64
LINX USB CHARGER PORT			
JBH2253		USB Charger Port DLX-USB02-A	£165
LINX USB CHARGER PORT BRACKET			
JBH2281		RH Mount	○
JBH2282		LH Mount	○
REMOTE STOP SWITCH			
JBH2258		Remote Stop Switch (ASL504)	(NC 2203, 2204) £1,050
LINX IN200 MODULE			
JBH2273		IN200 Module (Jack port for external power button)	£241
LINX INPUT MODULE			
UNI9252		IN500 Module (for DB9 Drive input)	(NC 2201, 2203, 2203) £785
FURTHER LINX COMPONENTS / ACCESSORIES			
JBH2300		To be enclosed only	○
JBX2301		To be connected and programmed (Special)	POA

BUDDY BUTTONS			
JBH1685		Red Buddy Button	£88
JBH1686		Yellow Buddy Button	£88
JBH1687		Green Buddy Button	£88
JBH1690		Blue Buddy Button	£88
EGG SWITCH			
JBH2254		Egg Switch, black	£88
PIKO BUTTONS			
JBH1711		Red Piko Button	£158
JBH1712		Yellow Piko Button	£158
JBH1713		Green Piko Button	£158
JBH1714		Orange Piko Button	£158
JBH1715		Blue Piko Button	£158
JBH1716		Black Piko Button	£158
BODYPOINT JOYSTICK KNOBS			
JBH2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A	£63
JBH2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A	£63
JBH2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£63
JBH2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£63
5:1 SWITCH ADAPTOR			
JBH1681		5:1 Switch Adapter (9-pin sub D to 5 jack ports)	£121
SHROUD COLOUR*			
JBH2871		Mid range - Black/Anthracite Shroud	○
JBH2872		Mid range - Silver Shroud	○
JBH2873		Mid range - Bengal Red Shroud	○
JBH2874		Mid range - Mercury Blue Shroud	○
JBH2875		Mid range - Midnight Purple Shroud	○
JBX0299		Premium - Individual Shroud Colour (includes coloured spreader bars on Modulite back) RAL number required Extended lead lime: + 15 days	£240
RIM COLOUR*			
JBH2882		Mid range - Black Rim Inserts	○
JBH2883		Mid range - Silver Rim Inserts	○
JBH2884		Mid range - Bengal Red Rim Inserts	○
JBH2885		Mid range - Mercury Blue Rim Inserts	○
JBH2886		Mid range - Midnight Purple Rim Inserts	○
JBX2887		Premium - Individual Rim Insert Colour RAL number required Extended lead lime: + 15 days	○
FORKS*			
JBH2631		Single sided fork	£226
JBH2632		Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆
FRONT WHEEL *			
JBH0835		9" (225 mm) x 2.80 (70 mm) Castor, pneumatic	○
JBH0836		9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	◆

CASTOR LOCKS			
JBH1490		Castor Lock	£197
DRIVE WHEEL*			
JBH0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
JBH0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	◆
C.T.C SUSPENSION*			
JBH2870		CTC Suspension Control - Traction - Comfort	◆
SUSPENSION SETTING*			
JBH1916		Suspension up to 160 kg User Weight For RX 40 HD only	◆
MOTOR TYPE*			
JBH1294		Indoor setting, reduced speed approx. 3 km/h	○
JBH1282		6 km/h	◆
JBH1908		10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 2203, 1437, 1254, 2837) £543
JBH1909		12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 2201, 2203, 1437, 1254, 2837) £653
BATTERY TYPE*			
JBH1259		60 Ah Gel	◆
JBH1260		73,5 Ah Gel	(NC 1450) £92
G TRAC			
JBH1215		G-Trac Module	£324
CHARGER*			
JBH1343		Opti Charge Charger 8A (IP44)	◆
JBH1313		Opti Charge Charger 10A (IP44)	£151
ANTI-TIPPERS*			
JBH1205		Anti-tipper	◆
LIGHTS			
JBH1437		LED Light system with Indicators 30 Lux headlights	(NC 2203, 1908, 1909) £479
JBH2837		LED Light system extra bright with indicators 60 Lux headlights	(NC 2203) £143
JBH1254		Without Light system	(NC 1909) ○
KERB RIDER			
JBH1244		Kerb Rider With park position feature	(NC 0556, 0589, 1294) £212
TRANSPORTATION KIT*			
JBH1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory) Includes a mounting bracket in the battery box Docking Station for car implementation has to be ordered separately with DAHL Engineering	£586
JBH1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory, Not compatible with Seat lifter	◆
JBH1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with combined Lift & Tilt module Adequate headrest is compulsory	(NC 1091, 1092, 1093) ○
EASY ACCESS FRONT TIE DOWNS			
JBH2894		Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down For transport worthy versions (1494/1401) only	£49

LAP BELT*			
JBH1329	<input type="checkbox"/>	Lap belt with Auto Lock Standard mount on backrest	○
JBH2113	<input type="checkbox"/>	Lap belt with Auto Lock Flexible mount in side rails	○
REAR VIEW MIRROR			
JBH1415	<input type="checkbox"/>	Rear view mirror RH	(NC 2217, 2285, 2213, 2265, 2275, 2283, 2287) £42
JBH1339	<input type="checkbox"/>	Rear view mirror LH	(NC 2218, 2286, 2214, 2266, 2276, 2284, 2288) £42
OTHER OPTIONS			
JBH2128	<input type="checkbox"/>	Bag Hooks on Modulite Seat Max load up to 6 kg = impact on max user weight	(NC 2227, 2228, 2233, 2240, 2257, 2291, 2292) £77
JBH1312	<input type="checkbox"/>	Cane holder RH	(NC 2227, 2259, 1221) £42
JBH1363	<input type="checkbox"/>	Cane holder LH	(NC 2228, 2260, 1222) £42
JBH1416	<input type="checkbox"/>	Fender for Castors	£88
JBH1417	<input type="checkbox"/>	Fender for Drive wheels	£169
JBH1404	<input type="checkbox"/>	Tray, Swing Away RH	(NC 1707, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2207, 2213, 2214, 2217) £167
JBH1405	<input type="checkbox"/>	Tray, Swing Away LH	(NC 1708, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2208, 2213, 2214, 2218) £167

TOTAL PRICE	
--------------------	--

Non Conforming Options:

JBX2171	Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm MUST BE SELECTED WITH HD SEATING JBX2102
JBX2197	Long Ergonomic Armpad (397 mm)
JBX1350	Neckrest depth adjustable
JBX1432	Headrest sideways adjustable
JBX1427	Headrest with cheek support
JBX2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware
JBX2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23) Complete with Elan Headrest Hardware
JBX2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard Complete with Elan Headrest Hardware
JBX1091	Manually settable Seat Tilt 0 - 15°
JBX1092	Powered Tilt (0° - 20°)
JBX1093	Powered Tilt (0° - 30°) with CoG shift
JBX1094	Seat lifter with powered tilt (300 mm, 30°) with CoG shift, 6 point Tie downs Must be selected
JBX0977	Powered Recliner 30°
JBH0556	Centre mount legrest, angle settable Length: 280 - 385 mm
JBH0589	LNx Power centre mount legrest, length compensating Length: 340 - 410 mm
JBH0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
JBH0503	Hanger Vari F pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with tools
JBH0533	Hanger Vari A pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with handwheel and calf pad
JBH0531	Hanger ADE pair Lower leg length: 290 - 460 mm Powered Angle adjustment: 80 - 0° with calf pad
JBH0504	Angle Adjustable Footplate pair
JBH0513	Angle and depth adjustable Footplate pair
JBH2201	Remote LiNX DLX-REM110 = Drive only 1 Drive function with Drive only power module
JBH2203	Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators
JBH2204	Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights
JBH2205	Rehab Remote LiNX, touch-screen user interface DLX-REM400 Configurable jack sockets
JBH2206	Rehab Remote LiNX without joystick for alternative access users DLX-REM500 Requires a user input/alternative joystick: 2231, 2232, 2229, 2230, 2247, 2233
JBH2270	Assignable Toggle Switch Module f. REM400 (DLX-SW400TGL)
JBH1207	Standard Bracket RH
JBH1208	Standard Bracket LH
JBH1707	Swing Away Bracket RH
JBH1708	Swing Away Bracket LH
JBH1709	Matrx Resolve Swing in/out bracket RH
JBH1710	Matrx Resolve Swing in/out bracket LH
JBH2213	Centre Mount, swing-away to RH, Nucleus
JBH2214	Centre Mount, swing-away to LH, Nucleus
JBH2227	Bracket RH, mounted for Attendant
JBH2228	Bracket LH, mounted for Attendant
JBH2220	Dual Control LiNX (DLX-REM050 as secondary remote)

JBH2219	Compact Dual Control LiNX (DLX-ACU200)
JBH2245	Invacare LiNX IDC Dual control
JBH2259	Bracket RH, Attendant
JBH2260	Bracket LH, Attendant
JBH1221	Bracket RH, detachable
JBH1222	Bracket LH, detachable
JBH2231	Proportional Compact Remote DLX-CR400
JBH2232	Proportional Compact Remote Low Force DLX-CR-400LF
JBH2207	RH Mount
JBH2208	LH Mount
JBH2265	Centre Mount, swing-away to RH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!
JBH2266	Centre Mount, swing-away to LH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!
JBH2230	Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M3
JBH2247	Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings
JBH2257	Four-Switch Proximity Array (ASL106)
JBH2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBH2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBH2233	MEC Proportional Joystick (Micro Extremity Control ASL130)
JBH2240	Proportional Heavy Duty Remote Joystick w mode switch (ASL133)
JBH2285	Tray Mount RH (100x150 mm)
JBH2286	Tray Mount LH (100x150 mm)
JBH2283	Centre Mount, swing-away to RH, Nucleus
JBH2284	Centre Mount, swing-away to LH, Nucleus
JBH2287	Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm
JBH2288	Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm
JBH2217	Bracket RH, mounted for user
JBH2218	Bracket LH, mounted for user
JBH2291	Bracket RH, mounted for attendant
JBH2292	Bracket LH, mounted for attendant
JBH2258	Remote Stop Switch (ASL504)
JBH2273	IN200 Module (Jack port for external power button)
JBH1294	Indoor setting, reduced speed approx. 3 km/h
JBH1908	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBH1909	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBH1437	LED Light system with Indicators 30 Lux headlights
JBH2837	LED Light system extra bright with indicators 60 Lux headlights
JBH1254	Without Light system
JBH1244	Kerb Rider With park position feature
JBH1450	Preparation for DAHL-Docking station (Adequate headrest is compulsory) Includes a mounting bracket in the battery box Docking Station for car implementation has to be ordered separately with DAHL Engineering
JBH1401	6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with combined Lift & Tilt module Adequate headrest is compulsory

JBH1415	Rear view mirror RH
JBH1339	Rear view mirror LH
JBH2128	Bag Hooks on Modulite Seat Max load up to 6 kg = impact on max user weight
JBH1312	Cane holder RH
JBH1363	Cane holder LH
JBH1416	Fender for Castors
JBH1404	Tray, Swing Away RH
JBH1405	Tray, Swing Away LH



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):