

Action[®] 2_{NG}

TRANSIT LITE



Lighter and more compact

Even better for travelling



Yes, you can.®

Take *Action*[®] and Go.

The new Action 2 NG Transit Lite is the newest chair from our award-winning Action range. Comfortable and reliable, Action chairs open a world of potential. Each chair in the range has a highly versatile core design with customisable options. The Action range is recognised simply as one you can count on.





Say hello to the new Action 2 NG Transit Lite

The new Action 2 NG Transit Lite is now lighter, more compact and super easy to handle. Its reduced frame size makes it even better for travelling and fitting into snug storage spaces.

Cross compatible with the range of Action 2 NG options, it's easy to adjust and maintain to suit your needs. Supporting you wherever you need to go.

Let's embrace new adventures.



1



Lighter and more compact.

We've reduced the weight and frame size, making it easier to use and push. Action 2 NG Transit Lite has our narrowest frame ever and a shorter wheelbase. Manoeuvre with less effort, even around tighter spaces. It's simply easier to get around.

*From 12.5 kg transport weight.

2



Extra comfort.

We've designed a new Y-shaped backrest for extra comfort. Smooth steering and spring brakes let you travel with confidence. A newly positioned parking brake is accessible while seated, helping to maintain levels of independence.

3



Folds in seconds.

The compact folding frame and backrest allows the Action 2 NG Transit Lite to fit snugly into any space. Removable legrests and armrests detach in seconds. Easy to fold down and set up, always ready for your next journey.

New braking system for attendant and user

Ergonomic position with 2 backrest heights

New lighter front and rear wheels

Seat handles designed for easier handling

Reversible padded armrest can be adjusted for comfort

Ease of transfer with removable swing-away legrests

Ideal for storage



What's new



Reduced weight and frame size



Easily accessible parking brake

Technical data

For more information about this product, including the user manual, please visit our website:

www.invacare.co.uk

	Seat width	Seat depth	Seat angle	Seat height	Backrest angle	Backrest height	Front frame angle
Action 2 NG Transit Lite	380 / 405 / 430 / 455 / 480 mm	400 / 430 mm	3°	485 / 510 mm	7°	430 / 510 mm	80°
	Length of footrests	Armrest height	Driving unit width	Total width, folded	Total length inc. Footrests	Total length without footrests	Total height
Action 2 NG Transit Lite	330 / 460 mm	195 - 285 mm	Seat width +100 mm	250 mm	995 mm	760 mm	975 mm
	Max. Safe slope	Radius	Max. User weight	Total weight	Mass of the heaviest part		
Action 2 NG Transit Lite	7°	870 mm	125 kg	13.2 kg	10.5 kg		

Frame colour



Anthracite black



Discover the
full Action range



Tel: +44 1656 776222
UK@invacare.com
www.invacare.co.uk
www.invacare.ie

Invacare Ltd
Unit 4 Pencoed Technology Park
Bridgend CF35 5AQ
UK



Yes, you can.®